

*English*



*Innovation in Spirometry Oximetry Telemedicine*


[www.spirometry.com](http://www.spirometry.com)  
[www.mirsmartone.com](http://www.mirsmartone.com)

# Global spirometer market: company share analysis

---

*MIR distributes its products in over 100 countries and it has been ranked among the top 5 spirometry major players in the world!*

*MIR is identified as **visionary leader** for the strength of product portfolio, distribution network across the globe and business strategy excellence. \**



**\* Sources (2019):**

*Global Spirometry Market, Analysis and Forecast to 2023*  
Secondary Research, Company Annual Reports, Press Releases



Item Code 910004

INTERNATIONAL  
PATENT

Volume Accuracy  
 $\pm 2,5\%$  or 50 ml



# FlowMIR™

## Disposable Turbine Flowmeter

MIR exclusive product  
International patent

*Each turbine, which includes a cardboard mouthpiece, has been individually factory tested with a computerized system*

*For each patient, after the Spirometry test, both turbine and mouthpiece are thrown away*



FlowMIR Dispensers (10 and 60 pcs)



Comfortable  
Packaging



Singularly Tested  
and Packed



Eliminates staff  
clean-up time



No Sterilisation



No Calibration



No Cross  
Contamination



No Antibacterial  
Filter



Always 100%  
accurate and hygienic



ATS/ERS Compliant  
for accuracy



Not affected by  
Vapour Condensation



Not affected by  
Ambient conditions



The Best Sensor  
for Spirometry

ADVANCED



## Spirolab™ Bluetooth®

7 inch Touchscreen  
All in one portable Desktop Spirometer  
with Oximetry option




Real-time tests with **wireless Bluetooth™** and **USB**

Fast and silent **built-in printer** with customizable  
printout format

Lightweight for easy mobility, **carrying case** included,  
long life **rechargeable battery**.

**GLI predicted sets** available in Stand-Alone mode  
and PC mode

**Pediatric Incentive** available on display



On screen paediatric incentive



Always included: professional PC software,  
with unlimited free update



ATS/ERS 2019,  
CE 0476, FDA 510(k)



Optional function available:  
Oximeter with adult or paediatric finger probe



Available with both **DISPOSABLE**  
or **REUSABLE** turbine flowmeter



**Spirometry test:** FVC, VC, IVC, MVV, PRE/  
POST Bronchodilator comparison with a  
wide range of selectable parameters



Direct connection to an external **USB printer**





ADVANCED



Detachable flowmeter head to facilitate Oximetry test during 6 Minute Walk Test

## Spirodoc™ Bluetooth® Touchscreen Spirometer and 3D Oximeter™: 6MWT, Sleep test, 24h SpO2 Holter

Real-time tests with **wireless Bluetooth™** and **USB**

Long life Rechargeable Battery

Unique feature: can be configured for both Doctor's Office or for Patient Home use



Fully featured stand alone Oximeter



Patient data entry



Always included: professional PC software,  
with unlimited free update



ATS/ERS 2019,  
CE 0476, FDA 510(k)



**Optional function available:**  
Oximeter with adult or paediatric finger probe



Available with both **DISPOSABLE**  
or **REUSABLE** turbine flowmeter



**Spirometry test:**  
FVC, VC, IVC, MVV, PRE/POST  
Bronchodilator comparison with a  
wide range of selectable parameters

ADVANCED



## Spirobank II™ Advanced

### Spirometer with Oximetry Option

Portable minilab for **pulmonologists and respiratory therapists**

Real-time tests with **wireless Bluetooth™** and **USB**

Rechargeable, **long-life battery**

**Traffic lights** for immediate test interpretation

**SpO2 and pulse rate** (optional) directly on the display (including **plethysmographic curve**)



**Always included: professional PC software, with unlimited free update**



**ATS/ERS 2019, CE 0476, FDA 510(k)**



**Available with both DISPOSABLE or REUSABLE turbine flowmeter**



**Spirometry test:** FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters



**Optional function available:**  
Oximeter with adult or paediatric finger probe

## Minispir™

### Computer-based Spirometer

Real time Flow/Volume and Volume/Time curves on your PC for a comprehensive Spirometry

Includes a wide range of selectable parameters

Power via USB, no internal memory, no display, no maintenance

Embedded **Temperature sensor** for BTPS conversion



**Always included: professional PC software, with unlimited free update**



**ATS/ERS 2019, CE 0476, FDA 510(k)**



**Available with both DISPOSABLE or REUSABLE turbine flowmeter**



**Spirometry test:** FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters

LIGHT



## Spirobank II™ Basic

Essential,  
easy to use Spirometer

*Ideal for family doctors,  
occupational medicine, screening*

*Real-time tests via USB*

*Rechargeable, long-life battery*

*Traffic lights for immediate test interpretation*

*Backl, high-resolution display*

*On screen results and curves preview*



Always included: professional PC software,  
with unlimited free update



ATS/ERS 2019,  
CE 0476, FDA 510(k)



Available with both DISPOSABLE  
or REUSABLE turbine flowmeter



**Spirometry test:** FVC, VC, IVC, IC, ERV,  
PRE/POST Bronchodilator comparison  
with a wide range of selectable parameters

## Minispir Light™

Simplified Computer-based  
Spirometer

*Real time COPD and Asthma screening  
have never been so intuitive and inexpensive!*

*Spirometry test interpretation*

*Power via USB, no internal memory, no display,  
no maintenance*

*Embedded Temperature sensor for BTPS  
conversion*



Always included: professional PC software,  
with unlimited free update



Available with DISPOSABLE turbine flowmeter



**Spirometry test:** FVC, VC, IVC, PRE/POST  
Bronchodilator comparison

## ADVANCED SMART



# Spirobank II<sup>TM</sup> Smart Bluetooth SMART

Complete Spirometer with Oximetry option for PC, iPad and Android Tablets


Real-time tests with **wireless Bluetooth<sup>TM</sup>** connection to Tablets, and **USB connection** to PC

Rechargeable, **long-life battery**

**SpO<sub>2</sub> and Pulse rate** (optional) directly on the display (including plethysmographic curve)


Highly innovative **Virtual Assistant** available on the App, to support the operator before, during and after the test.

**Share test results** via the App, in PDF or HL7 to a remote server.



Quanjer 2012-GLI

Predicted values with LLN and Z-score



Pediatric Incentive



Always included: professional PC software, with unlimited free update



Always included: professional App (MIR Spiro) for iPad and Android Tablets



ATS/ERS 2019, CE 0476, FDA 510(k)



Available with both **DISPOSABLE** or **REUSABLE** turbine flowmeter



**Spirometry test:** FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters



**Optional function available:** Oximeter with adult or paediatric finger probe



PERSONAL



## Spiritel™ Bluetooth®

4 in One for Home use:  
Touchscreen Spirometer, 3D Oximeter™,  
Triaxial Accelerometer and eDiary

3D Oximeter™ with triaxial accelerometer allows to correlate  
SpO2 Desaturation Events, patient Physical Activity  
and Body Position during data recording

The best solution for remote management in Asthma, COPD,  
Cystic Fibrosis and Lung Transplant

Ideal for Clinical Trials with customizable protocols

Long life Rechargeable Battery



Customizable eDiary



Always included: professional PC software,  
with unlimited free update



Optional function available:  
Oximeter with adult or paediatric finger probe



Spirometry test: FVC, FEV1, FEV1/FVC Ratio,  
FEV6, FEF, FET, FEF2575

PERSONAL SMART



# Spirobank Smart™



The simplest device  
for cost-effective spirometry


Real time test on Smartphone and  
via **Bluetooth Smart 4.0**

Suitable for **all ages from 5 to 93 years** and  
multi-ethnic groups (GLI predicted sets)


**Professional PDF report**, share with anyone you like  
via Email, whatsapp and other Apps

Test results and spirometry guidelines are **easy to understand** via the App.

**Software Development Kit** (iOS and Android)  
available for developers



Flow/Volume Curve



Professional PDF report



**Spirometry test:** PEF, FEV1, FVC, FEV1/FVC, FEF2575, FEV6, VEXT, DTPEF, FEF75, FEF25, FEF50



Available with both **DISPOSABLE**  
or **REUSABLE** single patient turbine flowmeter



Always included: personal App (MIR Spirobank)  
for iOS and Android Smartphones

ATS/ERS 2019, CE 0476, FDA 510(k)



PERSONAL SMART



# Spirobank Oxi™



MIR Spirobank

Now available with SpO2 (%) and Pulse (BPM)


Real time **plethysmographic curve** on Smartphone via Bluetooth Smart 4.0

Ready to Connect: designed for **Remote Patient Monitoring**, Homecare, Clinical Trials and 3<sup>rd</sup> party software integration.


**Professional PDF report**, share with anyone you like via Email, whatsapp and other Apps

Innovative reflectance **Pulse-Oximetry sensor** (Touch). Easy to use and accurate

**Software Development Kit** (iOS and Android) available for developers



Real-time plethysmographic curve



Flow/Volume Curve



**Spirometry test:** PEF, FEV1, FVC, FEV1/FVC, FEF275, FEF6, VEXT, DTPF, FEF75, FEF25, FEF50



**Oximetry test:** SpO2 (%), Pulse (BPM), Plethysmographic curve



Available with both **DISPOSABLE** or **REUSABLE** single patient turbine flowmeter



Always included: personal App (MIR Spirobank) for iOS and Android Smartphones



ATS/ESR 2019, CE0476, FDA Pending



## PERSONAL SMART



# Smart One™



SmartOne

## Peak Flow and FEV1 on your Smartphone and Tablet

Ideal in the **self management** of Asthma, COPD, Lung Transplant care, Cystic Fibrosis and for use in Clinical Trials  
Suitable for **all ages from 5 to 93 years** and multi-ethnic groups (GLI predicted sets)

**Easy-to-read** test results. Symptoms scoring and notes can be added to each test

Share results in pdf (WhatsApp, email etc.) and print directly.

**Software Development Kit** (iOS and Android) available for developers




**Spirometry test:** PEF, FEV1



Available with **REUSABLE** single patient turbine flowmeter



Always included: **personal App** (MIR Smart One) for Smartphone and Tablet



Traffic light health indicator for PEF, symptoms scoring and notes



Graphic trend of Peak Flow or FEV1



ATS/ESR 2019, CE0476, FDA 510(k) - O.T.C.





PERSONAL SMART



# Smart One OXI™



SmartOne

Now available with SpO2 (%) and Pulse (BPM)

Real time **plethysmographic curve** on Smartphone and Tablet via Bluetooth Smart 4.0

Innovative reflectance **Pulse-Oximetry sensor (Touch)**.  
Easy to use and accurate

Ready-to-Connect with **3<sup>rd</sup> party Apps** for Personal Care

**Share results in pdf** (WhatsApp, email etc.) and print directly.

**Software Development Kit (iOS and Android)** available for developers



**Spirometry test:** PEF, FEV1




**Oximetry test:** SpO2 (%), Pulse (BPM),  
Plethysmographic curve




Available with **REUSABLE**  
single patient turbine flowmeter



Always included: **personal App (MIR Smart One)**  
for Smartphone and Tablet



Real-time  
plethysmographic  
curve







Traffic light health  
indicator for PEF,  
symptoms scoring  
and notes










ATS/ESR 2019, CE0476, FDA Pending

## COMPARISON MATRIX - GENERAL FEATURES

	Name	SPIROLAB	SPIROBANK II BASIC	SPIROBANK II ADVANCED	SPIROBANK II SMART	MINISPIR LIGHT	MINISPIR	SPIRODOC
	Picture							
DEVICE	Dimensions	8.6×8.3×2" (220×210×51 mm)	6×2×1" (160×55×25 mm)	6×2×1" (160×55×25 mm)	6×2×1" (160×55×25 mm)	5.5×1.9×1" (142×49.7×26 mm)	5.5×1.9×1" (142×49.7×26 mm)	4.2×2.4×1.0" (101×48×16 mm)
	Display	7"	2.5"	2.5"	2.5"	Not Required	Not Required	3"
	Weight	51 Oz (1450 g)	5 Oz (140 g)	5 Oz (145 g)	5 Oz (145 g)	2.2 Oz (65 g)	2.2 Oz (65 g)	4.0 Oz (116 g)
	Carrying Case	Included	Canvas Pouch	Canvas Pouch	Canvas Pouch	No	Canvas Pouch	Canvas Pouch
	Device Type	Desktop	Pocket Size	Pocket Size	Pocket Size	PC Based	PC Based	Pocket Size
	Real Time PC	Yes	Yes	Yes	Yes	Yes	Yes	Yes
	Real Time Tablet	No	No	No	Yes	No	No	No
	Power Supply	Rechargeable battery	Rechargeable battery	Rechargeable battery	Rechargeable battery	USB Socket of PC	USB Socket of PC	Rechargeable battery
	Battery Charger	Included	Not Included	Not Included	Not Included	Not Required	Not Required	Not Included
	Thermal Printer	Internal	No	No	No	No	No	No
	Connectivity	USB, Bluetooth 2.0	USB	USB, Bluetooth 2.0	USB, Bluetooth 4.0	USB	USB	USB, Bluetooth 2.0
	USB to Printer	Yes	No	No	No	No	No	No
POST BD	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Memory Spiro	10,000 tests	10,000 tests	10,000 tests	10,000 tests	Uses PC Memory	Uses PC Memory	10,000 tests	
Pulse-Oximetry	Optional	No	Optional	Optional	No	No	Optional	
SENSOR	Sensor Type	Turbine	Turbine	Turbine	Turbine	Turbine	Turbine	Turbine
	Disposable Sensor	Yes	Yes	Yes	Yes	Yes	Yes	Yes
	Reusable Sensor	Yes	Yes	Yes	Yes	No	Yes	Yes
SOFTWARE	Software Name Always Included	WinspiroPRO	WinspiroPRO	WinspiroPRO	WinspiroPRO	Winspiro Light	WinspiroPRO	WinspiroPRO
	USB to PC	Yes	Yes	Yes	Yes	Yes	Yes	Yes
	Bluetooth to PC	Yes	No	Yes	No only to tablet	No	No	Yes
	Network Version	Optional	Optional	Optional	Optional	No	Optional	Optional

MIR products comply with the latest ATS/ERS standardization of spirometry (September 2019 update), Volume accuracy  $\pm 2.5\%$  or 50 ml








## COMPARISON MATRIX - SPIROMETRY PARAMETERS

	Name	SPIROLAB	SPIROBANK II BASIC	SPIROBANK II ADVANCED	SPIROBANK II SMART	MINISPIR LIGHT	MINISPIR	SPIRODOC
	Picture							
SPIROMETRY (STAND ALONE)	Parameters List Spiro	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrapolated Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/ FVC, FEV2, FEV2/ FVC, FEV3, FEV3/ FVC, FEV6, FEV1/ FEV6, FEV1/PEF, FEV1/FEV0.5, FVC, FIV1, FIV1/ FVC, PIF, FIF25, FIF50, FIF75, FEF50/FIF50, VC, IVC, IC, ERV, IRV, Rf, VE, VT, ti, tE, VT/tI, tE/TTOT, MVV (measured), MVV (calculated)	FVC, FEV1, FEV1%, PEF, FEF25-75, FET, Extrap. Volume, Lung Age, VC, IVC, IC, ERV	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrap. Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/ FVC, FEV2, FEV2/ FVC, FEV3, FEV3/ FVC, FEV6, FEV1/ FEV6, FIVC, FIV1, FIV1/FVC, PIF, FIF25, FIF50, FIF75, FEF50/ FIF50, VC, IVC, IC, ERV, IRV, Rf, VE, VT, ti, tE, VT/ tI, tE/TTOT, MVV (measured), MVV (calculated)	FVC, FEV1, FEV1%, FEV6, PEF, FEF25- 75%, FVC, Lung Age, VC, IVC	No	No	FVC, FEV1, FEV1/ FVC%, FEV3, FEV3/FVC%, FEV6, FEV1/FEV6%, PEF, FEF25%, FEF50%, FEF75%, FEF25%-75%, FET, Estimated Lung Age, Extr. Vol., FIVC, FIV1, FIV1/ FIVC%, PIF, VC, IVC, IC, ERV, FEV1/ VC%, VT, VE, Rf, ti, te, ti/t-tot, VT/ ti, MVV measured, MVV (calculated)
	E-diary	No	No	No	No	No	No	Yes
SPIROMETRY (PC BASED)	Parameters List Spiro	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrapolated Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/ FVC, FEV2, FEV2/ FVC, FEV3, FEV3/ FVC, FEV6, FEV1/ FEV6, FEV1/PEF, FEV1/FEV0.5, FVC, FIV1, FIV1/ FVC, PIF, FIF25, FIF50, FIF75,	FVC, FEV1, FEV1%, PEF, FEF25-75, FET, Extrap. Volume, Lung Age, VC, IVC, IC, ERV	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrap. Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/ FVC, FEV2, FEV2/ FVC, FEV3, FEV3/ FVC, FEV6, FEV1/ FEV6, FIVC, FIV1, FIV1/FVC, PIF, FIF25, FIF50, FIF75, FEF50/ FIF50, VC, IVC, IC, ERV, IRV, Rf, VE, VT, ti, tE, VT/tI, tE/TTOT, MVV (measured), MVV (calculated)	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrap. Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/ FVC, FEV2, FEV2/ FVC, FEV3, FEV3/ FVC, FEV6, FEV1/ FEV6, FIVC, FIV1, FIV1/FVC, PIF, FIF25, FIF50, FIF75, FEF50/ FIF50, VC, IVC, IC, ERV, IRV, Rf, VE, VT, ti, tE, VT/tI, tE/TTOT, MVV (measured), MVV (calculated)	FVC, FEV1, FEV1/ FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FVC, Lung Age, VC, IVC	FVC, FEV1, FEV1%, FEV6, PEF, FEF25-75, FIVC, Lung Age, VC, IVC	FVC, FEV1, FEV1/ FVC%, FEV3, FEV3/FVC%, FEV6, FEV1/FEV6%, PEF, FEF25%, FEF50%, FEF75%, FEF25%-75%, FET, Estimated Lung Age, Extr. Vol., FIVC, FIV1, FIV1/ FIVC%, PIF, VC, IVC, IC, ERV, FEV1/ VC%, VT, VE, Rf, ti, te, ti/t-tot, VT/ ti, MVV measured, MVV (calculated)
	Parameters List Spiro	No	No	No	FVC, FEV1, FEV1%, PEF, FEF25-75, FET, ELA, IVC, EVC	No	No	No
SPIROMETRY (ON THE APP)	Upgradable	Free Lifelong	Free Lifelong	Free Lifelong	Free Lifelong	Free Lifelong	Free Lifelong	Free Lifelong

MIR products comply with the latest ATS/ERS standardization of spirometry (September 2019 update), Volume accuracy  $\pm 2.5\%$  or 50 ml

## COMPARISON MATRIX - OXIMETRY PARAMETERS (OPTION)

### OXIMETRY

Name	SPIROLAB	SPIROBANK II BASIC	SPIROBANK II ADVANCED	SPIROBANK II SMART	MINISPIR LIGHT	MINISPIR	SPIRODOC
Picture							
Spot Test	Yes	No	Yes	Yes	No	No	Yes
Spot Test Param	SpO2 [Baseline, Min, Max, Mean], Pulse Rate [Baseline, Min, Max, Mean], T90, T89, T88, T5, Index [12s], SpO2 Events, Pulse Rate Events [Bradycardia, Tachycardia]	No	SpO2 [Baseline, Min, Max, Mean], Pulse Rate [Baseline, Min, Max, Mean], T90, T89, T88, T5, Index [12s], SpO2 Events, Pulse Rate Events [Bradycardia, Tachycardia]	SpO2 [Baseline, Min, Max, Mean], Pulse Rate [Baseline, Min, Max, Mean]	No	No	SpO2 [Baseline, Min, Max, Mean], Pulse Rate [Baseline, Min, Max, Mean], T90, T89, T88, T5, Index [12s], SpO2 Events, Pulse Rate Events [Bradycardia, Tachycardia]
Walk Test	No	No	-	-	No	No	Yes
Walk Test Param.	-	No	-	-	No	No	O2-Gap, Estimated distance, Distance walked, Predicted distance [Min, Standard], TA2%, [SpO2±2%], TA4% [ΔSpO2±4%], Time [Rest, Walking, Recovery], Desaturation Area/ Distance Optional data entry: Borg Dyspnea [Baseline, End, Change], Borg Fatigue [Baseline, End, Change], Arterial blood pressure [Systolic, Diastolic], Oxygen administered
Sleep Test	No	No	-	-	No	No	Yes
Sleep Test Param.	-	No	-	-	No	No	Body position, SpO2 Events, Desaturation index (ODI), Desaturation [Mean Value, Mean duration, Longest duration, Nadir Peak], ΔSpO2 [Min Drop, Max Drop], Total Pulse Variations, Pulse Rate Index, NOD89% [SpO2<89%; >5min], NOD4% [SpO2 Basale<4%; >5min], NOD90% [SpO2<90%; Nadir<86%; >5min]
Memory Pulse-Ox	900 Hours Pulse-Oximetry	No	900 Hours Pulse-Oximetry	900 Hours Pulse-Oximetry	No	No	900 Hours Pulse-Oximetry

MIR products comply with the latest ATS/ERS standardization of spirometry (September 2019 update), Volume accuracy  $\pm 2.5\%$  or 50 ml

## COMPARISON MATRIX - SPIROBANK SMART AND SMART ONE

	Name	SPIROBANK OXI	SPIROBANK SMART	SMART ONE	SMART ONE OXI
DEVICE	Picture				
	Dimensions	1.93x4.29x0.82" (49x109x21 mm)	1.93x4.29x0.82" (49x109x21 mm)	1.93x4.29x0.82" (49x109x21 mm)	1.93x4.29x0.82" (49x109x21 mm)
	Display	No	No	No	No
	Weight	2.14 Oz (60.7 g)	2.14 Oz (60.7 g)	2.14 Oz (60.7 g)	2.14 Oz (60.7 g)
	Carrying Case	N/A	N/A	N/A	
	Device Type	Smart	Smart	Personal Smart	Personal Smart
	Real Time PC	No	No	No	No
	Real Time Smart Phone	Yes	Yes	Yes	Yes
	Power Supply	2 batteries AAA 1.5V	2 batteries AAA 1.5V	2 batteries AAA 1.5V	2 batteries AAA 1.5V
	Battery Charger	N/A	N/A	N/A	N/A
	Thermal Printer	No	No	No	No
	Connectivity	Bluetooth 4.0	Bluetooth 4.0	Bluetooth 4.0	Bluetooth 4.0
SENSOR	Memory on Device	No	No	No	No
	Pulse- Oximetry	No	No	No	Yes
	Sensor Type	Turbine	Turbine	Turbine	Turbine
	Disposable	Yes	Yes	No	No
APP	Reusable	Yes	Yes	Yes	Yes
	Oxi	No	No	No	Touch Reflective
	App Name	MIR Spirobank	MIR Spirobank, iSpirometry	MIR Smart One	Mir Smart One
	Android	Yes	Yes	Yes	Yes
	iOs	Yes	Yes	Yes	Yes
	Memory Spiro on App	No	No	Yes	Yes
	Memory Oxi on App	No	No	No	Yes
	BT to Printer	Yes	Yes	Yes	Yes
PARAMETERS	Parameters Spiro	PEF, FEV1, FVC, FEV1/FVC, FEF2575, FEV6, VEXT, DTPEF, FEF75, FEF25, FEF50	PEF, FEV1, FVC, FEV1/FVC, FEF2575, FEV6, VEXT, DTPEF, FEF75, FEF25, FEF50	PEF, FEV1	PEF, FEV1
	Parameters Oxi	SpO2 (%), Pulse (BPM) Plethysmographic curve	No	No	SpO2 (%), Pulse (BPM) Plethysmographic curve
	Symptoms scoring	Yes	Yes	Yes	Yes

MIR products comply with the latest ATS/ERS standardization of spirometry (September 2019 update), Volume accuracy  $\pm 2.5\%$  or 50 ml

# Winspiro PRO™ & Winspiro PRO NET™

Complete and powerful software package for Specialists.  
Includes Telemedicine and Home Care Management

## Winspiro PRO™

Customizable PC software for Spirometry,  
Oximetry and Telemedicine projects

*FVC, VC, IVC, MVV, Bronchodilator Reversibility, Bronchial  
Challenge with FEV1-response curve with protocols for both  
Methacholine and Mannitol*

*Pædiatric Incentive for Spirometry*

*Home Care and Telemedicine management*

*Quanjer 2012-GLI Predicted values with LLN and Z-score*

*Patient Trend Charts for easy follow-up*

*Search engine for instant access to data*

*Specialized and customizable Printout*

*Wide choice of communication protocols ... and more*

*3D Oximeter® management, O2Gap Index,*

*Sleep Desaturation Analysis, 6 Minute Walk Test, Daytime*

*Physical Activity Report coupled with Desaturation Events*


## Winspiro PRO NET™

Network version of Winspiro PRO™

PC software for Spirometry,  
Oximetry and Telemedicine

*Complements the features of the client version,  
by providing the capability to share a single  
Database between different Network users*

Winspiro PRO™ is compatible with Windows: XP, VISTA, 7, 8 and 10.  
Winspiro PRO™ includes free updates and is provided with Spirolab™,  
Spirodac™, Spirobank II™, Spirolet™, Minispir™, Spirolab III™



Screenshot summary of all tests carried out



6 Minute Walk Test:  
Events and Physical Activity recording

Winspiro PRO NET™ includes free  
updates, and is provided on request

# MIR SDK



Go-To-Market-Toolkit  
Software Development Kit (SDK)  
For Windows and MacOS (PODS)

## MIR Software Development Kit

*Is the Software Development Kit to easily connect and handle MIR spirometers and oximeters into your Mac OS and Windows application. The developers will take advantage of the many libraries included in the SDK covering many different environments such as:*



### **USB and Bluetooth® communication**

*Real time values for spirometry and oximetry are provided to be plotted in a chart.  
USB Windows Certified Driver comes with the SDK.*



### **Spirometry Reference Values**

*A huge database of equations from many authors is available. For each parameter the MirPredicted library retrieves, not only the standard reference values but also the LLN, the ULN and the zScore values*



### **Spirometry Interpretation**

*Spirometry Interpretation library performs the Acceptability and Repeatability calculation (Quality Control Grade) and retrieves the Variability and the Automatic Interpretation for both the baseline spirometry and the post bronchodilator sessions, according to the ERS/ATS criteria.*



### **Oximetry Interpretation**

*Starting from the measured SpO2 and HeartRate "points" this library is able to perform both the generic and the specific analysis of the pulse oximetry, also using sophisticated algorithms for the detection of specific disease such as OSAS (Obstructive Sleep Apnea Syndrome), NOD (Nocturnal Oxygen Desaturation).*



### **Spirometry and Oximetry Charts**

*MirCharting library returns complete charts for any spirometry or oximetry layouts*



### **Spirometry and Oximetry Report in PDF**

*The SDK provides a professional and complete, ready to print report in PDF format.*

*Direct integration via USB with following MIR models:*

*Spirolab, Spirodoc, Minispir, Spirobank II Advanced, Spirobank II Basic*

Windows based

# Winspiro Light<sup>TM</sup>

Simple and intuitive for essential Spirometry test


Software for Real Time Spirometry on your PC

Simple and intuitive user interface with one main screen with all functions available

Ideal for primary care and occupational medicine

For each session produces a report with the 3 best tests including both Flow/Volume loop and Volume/Time curves that can be exported in .DOC or .PDF format

Wide range of selectable predicted values





# MIR Spiro™ android iOS

Simple and intuitive App for Spirometry and Oximetry test directly on your Tablets



## An expert beside you

Highly innovative **Virtual Assistant** supporting the operator with:

- ✓ recommendations for the correct execution of the test
- ✓ identification of the cause of errors during test execution
- ✓ automatic test interpretation

## Sharing is caring


Test can be printed, exported in pdf and sent by email.

Data can be transferred in HL7 format to a web service and are enabled for **EMR access**


**Desktop version** soon available for PC and Mac, including:

Unprecedented choice of **EMR Connectivity**: ready to connect with Epic, Athen Health, All Scripts, and GE Centricity among many others

Compliance with the **latest standards** in Spirometry, Occupational Health, Privacy and Security: ATS/ERS 2019, NIOSH-OSHA, GDPR, HIPAA



Didactic Guide



Virtual Assistant

Go-To-Market toolkit: **Software Development Kit (API)** available for App developers

# MIR Spirobank App **android** **iOS**



Real time Spirometry and Oximetry test on patient's Smartphone

## Main Features

- ✓ Spirometry and Oximetry test in one App
- ✓ Automatic connection via Bluetooth Smart BLE 4.0
- ✓ Suitable for all ages excluding newborns and infants
- ✓ GLI-2012 predicted sets for spirometry reference value




## Distinctive Features

- ✓ Acceptability Messages and Real time animation on Smartphone, to improve personal compliance during the test
- ✓ E-diary, symptoms and notes can be added for each test.
- ✓ Easy-to-read graphic trends for self-monitoring of spirometry and oximetry.
- ✓ Flow/Volume loop and Real time plethysmographic curve on screen

## Professional PDF Report

Share PDF report via whatsapp, e-mail and other applications. Pdf Report include: Flow/Volume and Flow/time curve (Best 3 trials), Quality Control Grade, Acceptable Trials, Variability of FEV1 and FVC, Pictograms

Go-To-Market toolkit: **Software Development Kit (API)** available for App developers



# Smart One App **android** **iOS**

Real time Spirometry and Oximetry test  
for Personal Care on Smartphone  
and Tablet




## Main Features

- ✓ Spirometry and Oximetry test in one App
- ✓ Automatic connection via Bluetooth Smart BLE 4.0
- ✓ Suitable for all ages excluding newborns and infants
- ✓ GLI-2012 predicted sets for spirometry reference value

## Distinctive Features

- ✓ Acceptability Messages and Real time animation on Smartphone and Tablet, to improve personal compliance during the test
- ✓ E-diary, symptoms and notes can be added for each test.
- ✓ Easy-to-read graphic trends for self-monitoring of spirometry and oximetry.
- ✓ Acceptability Messages, "Traffic Light" result interpretation for Peak Flow.



Go-To-Market toolkit: **Software Development Kit (API)** available  
for App developers

**Reusable Turbine**

Code 910002

**MIR Model** All Spiro

Reusable turbine with mesh.

**FlowMIR Turbine without mouthpiece**

Code 910001.

**MIR Model** All Spiro

Disposable turbine, compatible with Disposable Bacterial/Viral filter and Mouthpiece

**Paper adult disposable mouthpiece**

Code 910300

**MIR Model** All Spiro

Disposable, individually packed paper adult mouthpiece (CE marked).

**Reusable mouthpiece**

Code 910307

**MIR Model** All Spiro

Plastic adult reusable mouthpiece.

**Plastic paediatric reusable mouthpiece**

Code 910310

**MIR Model** All Spiro

Plastic paediatric reusable mouthpiece.

**Bacterial/Viral Filter**

Code 910302

**MIR Model** All Spiro

Bacterial/Viral Filter, disposable, packed in individual bag.

## Consumables / Accessories



### Paper roll for printer

Code 910350

**MIR Model** Spirolab

Standard thermal paper roll for printer.



### Plastic noseclip

Code 910320

**MIR Model** All Spiro

Plastic noseclip.



### Calibration syringe 3L

Code 919000

**MIR Model** All Spiro

3L calibration syringe.



### Oximetry finger sensor

Code 919024

**MIR Model** All devices with oximetry option

Reusable, ADULT.



### Oximetry finger sensor

Code 939005

**MIR Model** All devices with oximetry option

Reusable, PAEDIATRIC 5-45 kg.

# Spirometry for all!

- ✓ **Doctor and patient**

MIR offers groundbreaking solutions for both care givers and patients.

- ✓ **Clinical Trials & System Integrators**

MIR products can communicate with third parties systems and are fully customizable.

- ✓ **Mobile Health**

Patients can self monitor themselves through Apps on their Smart phones with the same accuracy of a clinical setting and keep track of their values.

## MIR Products Targets:



Spirolab



Spirobank II Advanced  
Spirobank II Smart



Spirodoc



Minispir

ADVANCED

*Pulmonology, Allergy, Cardiology, Geriatrics, Sleep Labs, Primary Care, Paediatrics, Occupational Medicine*



Spirobank II Basic



Minispir Light

LIGHT

*Primary Care, Paediatrics, Occupational Medicine, Sports Medicine, Internal Medicine, Anaesthesiology (pre-surgery), ERs*



Spirobank Smart



Spirobank Oxi



Smart One



Smart One Oxi

PERSONAL /SMART

*Pharma Companies (Clinical Trials), eHealthcare (Telemedicine), Home Care Monitoring, Neuromuscular Disease Management, Asthma Camps, Remote Areas Screening, Mobile Spirometry*



[www.spirometry.com](http://www.spirometry.com)  
[www.mirsmartone.com](http://www.mirsmartone.com)

## ITALY

### MIR Head Office

Via del Maggiolino, 125  
00155 Roma

**Tel.** +39 06 22 754 777  
**Fax** +39 06 22 754 785

[mir@spirometry.com](mailto:mir@spirometry.com)

## USA

### MIR USA, Inc.

5462 S. Westridge Drive  
New Berlin, WI 53151

**Phone** +1 (262) 565-6797  
**Fax** +1 (262) 364-2030

[mirusa@spirometry.com](mailto:mirusa@spirometry.com)

## FRANCE

### MIR Local Office

Jardin des Entreprises,  
290, Chemin de Saint Dionisy  
30980 LANGLADE (France)

**Phone** +33 (0)4 66 37 20 68  
**Fax** +33 (0)4 84 25 14 32

[mirfrance@spirometry.com](mailto:mirfrance@spirometry.com)

MIR reserves the right to modify the technical characteristics at any time