



Dimension

External diameter	Ø 33.3 mm
Internal diameter	Ø 30.5 mm
Length	41.5 mm
Thickness	2 mm
weight	20.29 g

Material

Dark part	Makrolon 2805
Deflector	LEXAN 121R
Alveolus	Synthetic sapphire
Blade	Audio mastering tape
Axis for blade	Corepoint® (Schneider)
Mesh	steel
Phthalates	free
Latex	free

Manufacturer

MIR s.r.l.

Classification

CE class	I
CND code	Z12150180
GMDN code	46906