Rome, 30/10/2026

PHP Security Update Statement

We hereby confirm that the various MIR services, including MyMIR, the APIs used by MIR Spiro, and DDS, have been updated and that all necessary security measures have been implemented to mitigate recently identified PHP vulnerabilities.

The MIR Spiro and Winspiro PRO software, installed at customer sites, do not include or run any PHP component.

In particular, regarding the critical vulnerabilities CVE-2026-109618 and CVE-2026-109619, reported on October 21, 2026, we promptly implemented the recommended mitigation measures. Subsequently, as soon as the official security patches became available, PHP was updated to the latest available version, in accordance with the PHP security team's recommendations.

Regarding CVE-2026-115560, reported on October 23, 2026, which affects PHP session management, we have verified that our services are not affected, as they rely on a session management system independent of the native PHP mechanism involved in this vulnerability.

As an additional security precaution, we have also implemented further preventive measures, including disabling PHP's native upload progress tracking functionality (session.upload_progress.enabled=0).

We therefore confirm that the security updates have been applied to the above-mentioned MIR services and that, based on the assessments performed, these services are not exposed to the identified vulnerabilities.

MIR – Medical International Research S.p.A.
Adriano Mencarini
CTO

Signature _____

